

Comprehensive insights into UV-Vis spectroscopy for the characterization and optical properties analysis of plasmonic metal nanoparticles

Giulia Conti¹, Lorenzo Moretti¹, Marco Rossi² and Elisa Ricci²

¹Department of Materials Science and Nanotechnology, University of Milan, Italy

²Department of Industrial Engineering, University of Naples Federico II, Italy

ABSTRACT

UV-Vis spectroscopy serves as a popular analytical method for characterizing plasmonic metal nanoparticles, thanks to its straightforwardness, high sensitivity, and capacity to offer real-time observations of their optical characteristics. The phenomenon of localized surface plasmon resonance (LSPR), which is a key feature of metal nanoparticles like gold and silver, is essential for grasping their distinctive optical behavior. This mini-review emphasizes the basic concepts of UV-Vis spectroscopy and its significant function in assessing nanoparticle size, shape, aggregation state, and stability. Important spectral characteristics, such as peak position, intensity, and full width at half maximum (FWHM), are examined concerning particle morphology and environmental conditions. Developments in the use of UV-Vis spectroscopy for tracking nanoparticle synthesis, stability, and interactions with biomolecules are also covered. Although UV-Vis spectroscopy is a powerful technique, it faces challenges like sensitivity to solvent effects and spectral overlap in mixed systems. Emerging trends, including the combination of UV-Vis with other characterization techniques and the use of machine learning for spectral interpretation, are discussed. This review highlights the flexibility and significance of UV-Vis spectroscopy in plasmonic nanoparticle research, providing insights relevant to both fundamental investigations and practical uses in areas such as sensing, catalysis, and biomedicine.

KEYWORDS

UV-Vis spectroscopy; Plasmonic metal nanoparticles; Optical properties analysis; Surface plasmon resonance; Nanoparticle characterization

ARTICLE HISTORY

Received 4 June 2024;
Revised 24 June 2024;
Accepted 2 July 2024

Introduction

Plasmonic metal nanoparticles, particularly those made of gold and silver, have gained considerable interest due to their distinctive optical characteristics that arise from localized surface plasmon resonance (LSPR) [1]. This effect, which results from the collective oscillation of conduction electrons when exposed to light, creates exceptional traits such as intense light absorption, scattering, and sensitivity to their environment. These qualities are fundamental to their extensive use in areas including biosensing, catalysis, photothermal therapy, and environmental monitoring [2].

UV-Vis spectroscopy is recognized as a robust, accessible, and economical method for analyzing plasmonic nanoparticles. Through the examination of peaks induced by LSPR in UV-Vis spectra, researchers can derive important information regarding particle size, shape, aggregation state, and surface modifications [3]. For example, alterations in the position, intensity, or width of these peaks offer direct insights into the physical and chemical attributes of the nanoparticles, making UV-Vis spectroscopy a crucial technique in the study of nanomaterials.

This presents a thorough investigation of UV-Vis spectroscopy as a method for the characterization and analysis of plasmonic metal nanoparticles. It covers how particle morphology affects spectral features, the function of UV-Vis in tracking nanoparticle stability, and its combination with other analytical techniques. Moreover, it emphasizes the challenges

and future prospects for enhancing the effectiveness of UV-Vis spectroscopy in nanoparticle research [4].

Fundamentals of UV-Vis Spectroscopy

UV-Vis spectroscopy is a commonly employed analytical technique for examining plasmonic metal nanoparticles by evaluating their absorption of ultraviolet and visible wavelengths. This approach is based on localized surface plasmon resonance (LSPR), where conduction electrons on the surface of the nanoparticle resonate with incoming light, producing characteristic absorption peaks. The location, strength, and breadth of these peaks are influenced by the size, shape, composition, and surrounding environment of the nanoparticles [5]. For instance, spherical gold nanoparticles typically display an LSPR peak near 520 nm, whereas anisotropic particles like rods or cubes may exhibit shifted or additional peaks. UV-Vis spectroscopy is also useful for tracking nanoparticle aggregation and stability, as aggregation results in peak broadening and red shifts. Its rapid, non-destructive analysis makes it particularly suitable for characterizing colloidal dispersions [6]. Given its user-friendly nature and versatility, UV-Vis spectroscopy is essential in the investigation and application of plasmonic nanomaterials.

UV-Vis Spectroscopy in Plasmonic Nanoparticle Characterization

UV-Vis spectroscopy is an extensively utilized and economical

*Correspondence: Dr. Giulia Conti, Department of Materials Science and Nanotechnology, University of Milan, Italy, e-mail: conti.g@yahoo.com

method for assessing plasmonic metal nanoparticles. It offers quick information about nanoparticle characteristics by leveraging the localized surface plasmon resonance (LSPR) effect, which affects their optical properties [7]. This method is crucial for evaluating the size, shape, aggregation, and stability of plasmonic nanoparticles.

Size and shape effects

The size and shape of plasmonic nanoparticles significantly affect their optical properties, especially their localized surface plasmon resonance (LSPR) peaks, which makes UV-Vis spectroscopy essential for their characterization [8] (Table 1):

- **Spherical nanoparticles:** Gold (520 nm) and silver (400 nm) nanoparticles, which are spherical in shape, display a single, sharp LSPR peak. As the size of the particles increases, the peak shifts to longer wavelengths (redshifts) and becomes broader.
- **Nonspherical nanoparticles:** Anisotropic shapes, such as rods, cubes, or triangles, present multiple LSPR peaks. For instance, gold nanorods have distinct peaks corresponding to the transverse (short axis) and longitudinal (long axis) plasmon resonance.

Aggregation and stability analysis

UV-Vis spectroscopy serves as an effective method for assessing the stability and aggregation status of plasmonic nanoparticles.

- **Effects of aggregation:** The aggregation of nanoparticles leads to a shift in their LSPR peaks towards longer wavelengths (redshift) and causes considerable broadening. This redshift arises from the plasmonic interaction between adjacent particles [9,10].
- **Evaluating stability:** By capturing spectra over time, one can assess the stability of nanoparticles within various media (such as solvents and buffers). This stability is vital for applications in fields like biomedicine and sensing, where nanoparticles must maintain their characteristics across different environments [11].

Advanced Applications

UV-Vis spectroscopy is essential in numerous applications involving plasmonic metal nanoparticles. In biosensing, it facilitates real-time and label-free identification of biomolecules by observing shifts in the localized surface plasmon resonance (LSPR) peak [12]. In environmental monitoring, it serves to identify harmful substances, while in catalysis, it assists in examining the stability and surface interactions of nanoparticles. The biomedical sector employs UV-Vis spectroscopy for the characterization of nanoparticles used in techniques like photothermal therapy and drug delivery [13,14]. Integrating UV-Vis with other methods, such as Raman

spectroscopy and electron microscopy, enhances the study of nanomaterials. Progress in machine learning and AI is also advancing the automated analysis of UV-Vis spectra, allowing for more accurate large-scale applications [15].

Challenges and Limitations

While UV-Vis spectroscopy is an effective method for analyzing plasmonic metal nanoparticles, it presents several challenges:

- **Sensitivity to environmental influences:** The spectral characteristics are affected by variables such as the polarity of the solvent, the refractive index, and temperature, which complicates data analysis [16].
- **Aggregation and size variation:** The presence of aggregated or polydisperse nanoparticles can lead to overlapping spectra, complicating the analysis of size and shape.
- **Interference from other components:** Impurities or other substances within complex samples can mask absorption peaks, diminishing measurement accuracy.
- **Limited information on structure:** UV-Vis spectroscopy primarily examines optical properties and does not provide direct insights into structural details, necessitating the use of additional techniques like TEM or XRD [17].
- **Challenges in quantitative analysis:** Scattering effects can obstruct accurate readings, particularly for larger or irregularly shaped nanoparticles.

To overcome these challenges, careful sample preparation, optimal conditions, and the integration of other analytical techniques are essential.

Future Perspectives

The prospects of UV-Vis spectroscopy in analyzing plasmonic metal nanoparticles center on enhancing sensitivity and adaptability. Innovations in instruments, such as advanced spectrometers and real-time monitoring technologies, will facilitate more precise and quicker identification of changes in nanoparticle properties [18,19]. Combining UV-Vis with methodologies like Raman spectroscopy, electron microscopy, and mass spectrometry will provide additional insights into the structure and chemistry of nanoparticles. The employment of machine learning for data interpretation holds the potential for quicker, automated analysis, tackling issues related to polydispersity and aggregation [20]. Furthermore, extending UV-Vis to investigate nanoparticles in dynamic settings, such as biological environments, will enrich our comprehension of their behavior in realistic scenarios, thereby improving applications in biomedicine, environmental surveillance, and energy conversion. These advancements will ensure that UV-Vis spectroscopy remains a leading technique in nanomaterials research [21].

Table 1. Relationship between Nanoparticle Morphology and UV-Vis Spectra

Nanoparticle Shape	Peak Position (LSPR)	Peak Broadening	Examples
Spherical	~520 nm (gold)	Narrow	Gold nanoparticles
Rod-shaped	~700–900 nm (long axis)	Broader	Gold nanorods
Cube-shaped	~550–600 nm	Moderate	Silver nanocubes

Conclusion

UV-Vis spectroscopy is a fundamental technique for characterizing plasmonic metal nanoparticles, providing essential information about their optical characteristics and behavior. The method's capacity to identify and assess the localized surface plasmon resonance (LSPR) delivers vital data regarding particle size, shape, aggregation, and interactions with the surrounding medium. It is crucial to comprehend the spectral characteristics, such as peak location, intensity, and broadening, to enhance nanoparticle performance across various applications, including sensing, catalysis, and biomedical treatments.

Despite its strengths, UV-Vis spectroscopy encounters constraints related to environmental conditions, sample variability, and difficulties in deciphering complex spectra. However, integrating UV-Vis with additional analytical methods, like electron microscopy or Raman spectroscopy, improves the depth and precision of nanoparticle characterization.

In the future, innovations in UV-Vis spectroscopy, including combined technologies and data-driven methodologies, are expected to advance the field of nanoparticle research, allowing for more accurate and effective analyses. This mini-review highlights the continued relevance of UV-Vis spectroscopy as a critical tool in the realm of plasmonic nanomaterials, assisting in the advancement of next-generation nanodevices and their applications.

Disclosure Statement

No potential conflict of interest was reported by the authors.

References

1. Hang Y, Wang A, Wu N. Plasmonic silver and gold nanoparticles: shape-and structure-modulated plasmonic functionality for point-of-caring sensing, bio-imaging and medical therapy. *Chem Soc Rev.* 2024;53:2932-2971. <https://doi.org/10.1039/D3CS00793F>
2. Jana J, Ganguly M, Pal T. Enlightening surface plasmon resonance effect of metal nanoparticles for practical spectroscopic application. *RSC Adv.* 2016;6(89):86174-86211. <https://doi.org/10.1039/C6RA14173K>
3. Pashkov DM, Guda AA, Kirichkov MV, Guda SA, Martini A, Soldatov SA, et al. Quantitative analysis of the UV-vis spectra for gold nanoparticles powered by supervised machine learning. *J Phys Chem C.* 2021;125(16):8656-8666. <https://doi.org/10.1021/ACS.jpcc.0c10680>
4. Ray TR, Lettiere B, de Rutte J, Pennathur S. Quantitative characterization of the colloidal stability of metallic nanoparticles using UV-Vis absorbance spectroscopy. *Langmuir.* 2015;31(12):3577-3586. <https://doi.org/10.1021/la504511j>
5. Mcoyi MP, Mpofo KT, Sekhwama M, Mthunzi-Kufa P. Developments in localized surface plasmon resonance. *Plasmonics.* 2024;1-40. <https://doi.org/10.1007/s11468-024-02620-x>
6. Alin J, Rubino M, Auras R. Effect of the solvent on the size of clay nanoparticles in solution as determined using an ultraviolet-visible (UV-Vis) spectroscopy methodology. *Appl Spectrosc.* 2015;69(6):671-678. <https://doi.org/10.1366/14-07704>
7. Kim J, Lee SK, Schellingerhout D, Nahrendorf M, Kim K, Kim J, et al. Spectroscopic Assessment of Gold Nanoparticle Biodistribution Using Surface Plasmon Resonance Phenomena. *ACS Biomater Sci Eng.* 2019;5(12):6389-6394. <https://doi.org/10.1021/acsbomaterials.9b01079>
8. Olson J, Dominguez-Medina S, Hoggard A, Wang LY, Chang WS, Link S. Optical characterization of single plasmonic nanoparticles. *Chem Soc Rev.* 2015;44(1):40-57. <https://doi.org/10.1039/c4cs00131a>
9. Li Z, Wang W, Yin Y. Colloidal assembly and active tuning of coupled plasmonic nanospheres. *Trends Chem.* 2020;2(7):593-608. <https://doi.org/10.1016/j.trechm.2020.03.008>
10. Kar D, Sen P, Srinivas V, Kasiviswanathan S. Effect of interparticle interaction on the plasmon resonance of silver nanoparticles. *AIP Conf Proc.* 2018;1942(1). <https://doi.org/10.1063/1.5028740>
11. Carvalho PM, Felício MR, Santos NC, Gonçalves S, Domingues MM. Application of light scattering techniques to nanoparticle characterization and development. *Front Chem.* 2018;6:237. <https://doi.org/10.3389/fchem.2018.00237>
12. Begum R, Farooqi ZH, Naseem K, Ali F, Batool M, Xiao J, et al. Applications of UV/Vis spectroscopy in characterization and catalytic activity of noble metal nanoparticles fabricated in responsive polymer microgels: A review. *Crit Rev Anal Chem.* 2018;48(6):503-516. <https://doi.org/10.1080/10408347.2018.1451299>
13. Borowska M, Jankowski K. Basic and advanced spectrometric methods for complete nanoparticles characterization in bio/eco systems: current status and future prospects. *Anal Bioanal Chem.* 2023;415(18):4023-4038. <https://doi.org/10.1007/s00216-023-04641-7>
14. Pandit J, Alam MS, Ansari JR, Singhal M, Gupta N, Waziri A, et al. Multifaceted applications of nanoparticles in biological science. *In Nanomaterials in the battle against pathogens and disease vectors.* 2022;17-50. <https://doi.org/10.1201/9781003126256-2>
15. Botifoll M, Pinto-Huguet I, Arbiol J. Machine learning in electron microscopy for advanced nanocharacterization: current developments, available tools and future outlook. *Nanoscale Horiz.* 2022;7(12):1427-1477. <https://doi.org/10.1039/d2nh00377e>
16. Monteiro AD, Feital TD, Pinto JC. Statistical aspects of near-infrared spectroscopy for the characterization of errors and model building. *Appl Spectrosc.* 2017;71(7):1665-1676. <https://doi.org/10.1177/0003702817704587>
17. Cole JM, Gosztoła DJ, Sylvester SO, Wang SG, Chen YS. Assigning optical absorption transitions with light-induced crystal structures: case study of a single-crystal nano-optomechanical transducer. *J Phys Chem C.* 2021;125(28):15711-15723. <https://doi.org/10.1021/ACS.jpcc.1c04425>
18. Krajczewski J, Kołataj K, Kudelski A. Plasmonic nanoparticles in chemical analysis. *RSC Adv.* 2017;7(28):17559-17576. <https://doi.org/10.1039/C7RA01034F>
19. Mohammadparast F, Tirumala RT, Ramakrishnan SB, Dadgar AP, Andiappan M. Operando UV-Vis spectroscopy as potential in-line PAT system for size determination of functioning metal nanocatalysts. *Chem Eng Sci.* 2020;225:115821. <https://doi.org/10.1016/j.ces.2020.115821>
20. Lu S, Jayaraman A. Machine learning for analyses and automation of structural characterization of polymer materials. *Prog Polym Sci.* 2024;101828. <https://doi.org/10.1016/j.progpolymsci.2024.101828>
21. Zamora-Perez P, Tsoutsis D, Xu R, Rivera_Gil P. Hyperspectral-enhanced dark field microscopy for single and collective nanoparticle characterization in biological environments. *Mater.* 2018;11(2):243. <https://doi.org/10.3390/ma11020243>